



AESCULAP CASE CART, WASH BIN, STORAGE BIN

INTRODUCTION

This Aesculap case cart, wash bin and storage bins are intended to hold reusable instruments, single use surgical packages and accessories, as assemblies during washing, storage and transport to/from the Operating Room. Unless explicitly stated, these products are NOT for transport outside a healthcare facility.

CARE AND HANDLING

Do not use abrasive cleaners, metal brushes, or abrasive cleaning pads. Use of abrasive products can cause permanent damage to surfaces. Use of abrasive cleanroom pads will result in warranty voidance.

INSPECTION PRIOR TO USE

Product should be free from:

- Noticeable cracking, corrosion, pitting and damage
- Any impairment to which cleaning methods not adequately made

Material conditions are not met, do not use, please report to supplier

- All components such as wheels and casters should move with appropriate ease - if not tolerate with suitable roller report

DISINFECTANTS AND DETERGENTS

Use only mild neutral detergents/low, neutral pH (7) detergents to clean effectively without causing damage to the surface. If in doubt, contact the detergent manufacturer or supplier to determine suitability for use in cleaning the Aesculap product and material.

Disinfecting Instructions

1. Remove all visible bio burden and debris

Manual Cleaning

1. Use a soft sponge and a mild detergent and clean the product per supplier's advice
2. Rinse thoroughly to remove all detergent residue
3. Thoroughly dry all components with a soft dry cloth

Always wear proper protective personal attire when cleaning

Mechanical Cleaning

1. Soak items up, if applicable, place the product in the washer, let drain if possible to avoid water collection
2. Soak reusable status and other items, if applicable, to minimize damage to washer and product
3. Thoroughly dry (either with a soft, dry cloth or using the components before use)

If you have further questions regarding cleaning practices, contact your local Aesculap Sales Representative or the Customer Service Dept. (toll free) at 1-800-352-5000

ASSEMBLY FROM KIT

Reactive Resin such as light curing composites, wrapped instrument sets, post practices only. No active Resin are present on or in the product, for applicable use. Please refer to 010412-01 for additional information on loading if needed.

Aesculap Inc.
3000 Lakeside Blvd. Ste. 300
Mettuchen, NJ 07043
Customer Service: 1-800-351-5500

AESCLAP

Aesculap Inc.
3000 Lakeside Avenue
Mettuchen, NJ 07043

010412-01

010412-01 Rev. 02

Effective Date: 03/15/12